

Today's Date: _____ Right Left Bilateral

PT Last Name _____

PT First Name _____

DOB _____ Weight _____ Sex M / F

DX _____

Company Name _____

Contact _____

Phone Number _____

Fax Number _____

Billing Address _____

City _____ ST _____ Zip _____

Shipping Address _____

City _____ ST _____ Zip _____

Email _____

How do you prefer to receive your order confirmation?

- FAX EMAIL

Pattern Options (Circle one)

- | | | | | |
|------------|---------------|--------------|----------|------------|
| White | BusyBees | Skateboard | Black | Royal Blue |
| Yellow | Butterflies | Sports | Lt .Blue | Pink |
| Blue | Construction | Traffic Soup | Yellow | Purple |
| Pink | Dinosaurs | Brown Camo | Red | White |
| BabySports | FarmFriends | Pink Camo | Green | |
| Beach | HappyDaisies | Pink Swirls | | |
| Blossom | Jungle Juggle | | | |
| Bugs | Ocean Life | | | |

Strap Color (Circle one)

- Chafe Option** (Circle one)
- Plastic (Standard)
- Dacron (Optional, no charge)

SureStep Children's Footwear



Please ship shoes with this order




Quantity (prs.) White/Pink White/Blue


Shoes are available in infant size 3 through toddler size 12.
(Measurement #6 = 2^{7/8}" - 5^{3/8}")

For volume shoe orders, please call **877-462-0711**
or visit www.surestepchildrensfootwear.com for details.


12 Circumference at distal fibula head




Smallest supra-malleolar circumference




Diagonal Circumference at apex of calcaneus



Midfoot circumference at apex of arch



11 Posterior heel to finished lateral length



LATERAL


MEDIAL

7

6



Posterior calcaneus to apex of 5th met-head




Posterior calcaneus to apex of 1st met-head

○ = Circumference

△ = Width

□ = Distance

* Width at met-heads (simulated weight bearing)



* If 2^{3/4}" or greater, a cast is required

10

Posterior heel to finished medial length

* Shoes are available for feet measuring 2^{7/8}" - 5^{3/8}"



Width at apex of malleoli

LATERAL

MEDIAL

Height to apex of lateral malleolus

Height to apex of medial malleolus

9

8

13 Desired finished height of AFO



SureStep AFO 2011

FAX ORDERS TO: 866-700-7837

Choose One: All Appropriate Options

advanced (solid ankle)

PullOver (hinged)

Dorsiflexion assist

2Stage

Hinged

Solid Ankle

Plantar flexion stop

Footplate (choose one)

- Full Length-standard
*10 & 11 required
- SureStep Trimlines

Other Options

- Metatarsal Arch Support
- Toe Plateau (to install at fitting)

Special Modifications:

Shipping Options (Circle one)

UPS Ground- No Charge

2nd day by Noon- \$17.00

Next day by 10:30a.m.- \$35.00

24 HOUR RUSH! ADD \$50 FEE
(SHIPPING CHARGES APPLY SEPARATELY)

NEED BY DATE: _____

SEND CASTS TO:

17530 Dugdale Drive
Phone: 1-877-462-0711

South Bend, IN 46635
Fax: 1-866-700-7837