

Today's Date: _____ Right Left Bilateral

PT Last Name _____

PT First Name _____

DOB _____ Weight _____ Sex M / F

DX _____

Company Name _____

Contact _____

Phone Number _____

Fax Number _____

Billing Address _____

City _____ ST _____ Zip _____

Shipping Address _____

City _____ ST _____ Zip _____

Email _____

PO _____

How do you prefer to receive your order confirmation?

- FAX EMAIL

Pattern Options (Circle one)

- | | | |
|------------|---------------|--------------|
| White | BusyBees | Skateboard |
| Yellow | Butterflies | Sports |
| Blue | Construction | Traffic Soup |
| Pink | Dinosaurs | Brown Camo |
| BabySports | FarmFriends | Pink Camo |
| Beach | HappyDaisies | Pink Swirls |
| Blossom | Jungle Juggle | |
| Bugs | Ocean Life | |

Strap Color (Circle one)

- | | |
|----------|------------|
| Black | Royal Blue |
| Lt .Blue | Pink |
| Yellow | Purple |
| Red | White |
| Green | |

Chafe Option (Circle one)

- Plastic (Standard)
- Dacron (Optional, no charge)

SureStep Children's Footwear




Please ship shoes with this order

Quantity (prs.)	White/Pink	White/Blue
_____	_____	_____

Shoes are available in infant size 3 through toddler size 12.
(Measurement #6 = 27/8" - 53/8")


For volume shoe orders, please call **877-462-0711**.

12 Circumference at distal fibula head




Smallest supra-malleolar circumference

1




2

Diagonal Circumference at apex of calcaneus



3

Midfoot circumference at apex of arch



SureStep AFO 2011

--FAX ORDERS TO: 866-700-7837--

Choose One: All Appropriate Options

advanced (solid ankle)

PulOver (hinged)

Dorsiflexion assist

Indy 2Stage

Hinged

Solid Ankle

Plantar flexion stop

Footplate (choose one)

Full Length-standard

*10 & 11 required

SureStep Trimlines

Other Options

Metatarsal Arch Support

Toe Plateau (to install at fitting)

4 day fab time

*Cast + measurements required

11 Posterior heel to finished lateral length



LATERAL

7

MEDIAL

6

*If 2 3/4" or greater, a cast is required


10

Posterior heel to finished medial length

*Shoes are available for feet measuring 27/8" - 53/8"



Posterior calcaneus to apex of 5th met-head



Posterior calcaneus to apex of 1st met-head



○ = Circumference

△ = Width

□ = Distance

4

Width at apex of malleoli

LATERAL

MEDIAL

Height to apex of lateral malleolus


9

Height to apex of medial malleolus

8

13

Desired finished height of AFO



Special Modifications:

Shipping Options (Circle one)

UPS Ground- No Charge

2nd day by Noon- \$17.00

Next day by 10:30a.m.- \$35.00

24 HOUR RUSH! ADD \$50 FEE

(SHIPPING CHARGES APPLY SEPARATELY)

NEED BY DATE: _____

SEND CASTS TO:

17530 Dugdale Drive
Phone: 1-877-462-0711

South Bend, IN 46635
Fax: 1-866-700-7837